

Product Datasheet

HIV 1 p24 (1D1) monoclonal antibody (orb1726146)

Description	HIV 1 p24 (1D1) monoclonal antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant full-length HIV-1 p24 protein.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P04591
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.