

Product Datasheet

Influenza A Hemagglutinin H5 (2C10) monoclonal antibody (orb1726148)

Description	Influenza A Hemagglutinin H5 (2C10) monoclonal antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified avian influenza virus type A (H5N1).
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q77XR3
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.