

## Product Datasheet

# Anti-Rabbit/Mouse EZ-HRP; Polymer Detection (orb1726438)

<b>Description</b>	Anti-Rabbit/Mouse EZ-HRP; Polymer Detection
<b>Reactivity</b>	Mouse, Rabbit
<b>Tested Applications</b>	IHC
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	4°C for up to 1 year
<b>Note</b>	For research use only
<b>Dilution Range</b>	0.458333333
<b>Expiration Date</b>	12 months from date of receipt.