

Product Datasheet

ErbB2(Phospho-Tyr1221 / Tyr1222) Rabbit mAb Antibody (orb1726549)

Description	ErbB2(Phospho-Tyr1221 / Tyr1222) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.