

## Product Datasheet

# CaMKII alpha (Phospho-Thr286) Rabbit mAb Antibody (orb1726552)

<b>Description</b>	CaMKII alpha (Phospho-Thr286) Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Recombinant Rabbit
<b>Uniprot ID</b>	<b>P11798, P11275, Q9UQM7</b>
<b>Expiration Date</b>	12 months from date of receipt.