


Product Datasheet

MYH7B Rabbit mAb Antibody (orb1726833)

Description	MYH7B Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Immunogen sequence is 1750-1983 amino acids of human MYH7B
Target	MYH7B Antibody detects endogenous levels of total MYH7B
Form/Appearance	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Purity	Purified by affinity-chromatography
Source	Recombinant Rabbit
MW	221 kDa
Uniprot ID	A7E2Y1
Expiration Date	12 months from date of receipt.



Western blot analysis of MYH7B expression...