

Product Datasheet

CNTN2 / TAG1 Rabbit mAb Antibody (orb1727187)

Description	CNTN2 / TAG1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	Q02246
Expiration Date	12 months from date of receipt.