

Product Datasheet

RUNX1/RUNX2/RUNX3 Rabbit mAb Antibody (orb1727594)

Description	RUNX1/RUNX2/RUNX3 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	Q13761 , Q13950 , Q01196
Expiration Date	12 months from date of receipt.