

## Product Datasheet

# PBK/TOPK (Phospho-Thr9) Rabbit mAb Antibody (orb1727670)

<b>Description</b>	PBK/TOPK (Phospho-Thr9) Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Storage</b>	Store at -20°C. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Recombinant Rabbit
<b>Uniprot ID</b>	<a href="#">Q96KB5</a>
<b>Expiration Date</b>	12 months from date of receipt.