

Product Datasheet

PDHA1 (Phospho-Ser293) Rabbit mAb Antibody (orb1727692)

Description	PDHA1 (Phospho-Ser293) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P08559
Expiration Date	12 months from date of receipt.