

## Product Datasheet

# IKB alpha (Phospho-Ser36) Rabbit mAb Antibody (orb1727693)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | IKB alpha (Phospho-Ser36)<br>Rabbit mAb Antibody |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse                                     |
| <b>Conjugation</b>         | Unconjugated                                     |
| <b>Tested Applications</b> | FC, WB   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze/thaw cycles.        |
| <b>Note</b>                | For research use only                            |
| <b>Clonality</b>           | Recombinant                                      |
| <b>Source</b>              | Recombinant Rabbit                               |
| <b>Uniprot ID</b>          | <b>P25963</b>                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                  |