

Product Datasheet

JAK2 (Phospho-Tyr1007/Tyr1008) Rabbit mAb Antibody (orb1727714)

Description	JAK2 (Phospho-Tyr1007/Tyr1008) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	O60674
Expiration Date	12 months from date of receipt.