

Product Datasheet

B Raf (Phospho-Thr401) Rabbit mAb Antibody (orb1727722)

| | |
|----------------------------|---|
| Description | B Raf (Phospho-Thr401) Rabbit mAb Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IP, WB |
| Storage | Store at -20°C. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Recombinant |
| Source | Recombinant Rabbit |
| Uniprot ID | P15056 |
| Expiration Date | 12 months from date of receipt. |