

Product Datasheet

Rac1/Cdc42 (Phospho-Ser71) Rabbit mAb Antibody (orb1727724)

Description	Rac1/Cdc42 (Phospho-Ser71) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P63000, P60953
Expiration Date	12 months from date of receipt.