

## Product Datasheet

# Cdc6 (Phospho-Ser54) Rabbit mAb Antibody (orb1727728)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Cdc6 (Phospho-Ser54)<br>Rabbit mAb Antibody |
| <b>Species/Host</b>        | Rabbit                                      |
| <b>Reactivity</b>          | Human                                       |
| <b>Conjugation</b>         | Unconjugated                                |
| <b>Tested Applications</b> | ICC, WB                                     |
| <b>Storage</b>             | Store at -20°C. Avoid freeze/thaw cycles.   |
| <b>Note</b>                | For research use only                       |
| <b>Clonality</b>           | Recombinant                                 |
| <b>Source</b>              | Recombinant Rabbit                          |
| <b>Uniprot ID</b>          | <a href="#">Q99741</a>                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.             |