

Product Datasheet

AKT1 (Phospho-Ser124) Rabbit mAb Antibody (orb1727729)

Description	AKT1 (Phospho-Ser124) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P31749
Expiration Date	12 months from date of receipt.