

Product Datasheet

EGFR (Phospho-Tyr1173) Rabbit mAb Antibody (orb1727773)

Description	EGFR (Phospho-Tyr1173) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.