

Product Datasheet

p53 (Phospho-Thr55) Rabbit mAb Antibody (orb1727790)

Description	p53 (Phospho-Thr55) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P04637
Expiration Date	12 months from date of receipt.