

## Product Datasheet

# Glycogen synthase 1 (Phospho-Ser641) Rabbit mAb Antibody (orb1727795)

<b>Description</b>	Glycogen synthase 1 (Phospho-Ser641) Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, IP, WB
<b>Storage</b>	Store at -20°C. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Recombinant Rabbit
<b>Uniprot ID</b>	<b>P13807</b>
<b>Expiration Date</b>	12 months from date of receipt.