

Product Datasheet

CDK1/2/3 (Phospho-Thr14) Rabbit mAb Antibody (orb1727800)

Description	CDK1/2/3 (Phospho-Thr14) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P24941, P06493, Q00526
Expiration Date	12 months from date of receipt.