

Product Datasheet

alpha Synuclein (Phospho-Ser129) Rabbit mAb Antibody (orb1727812)

Description	alpha Synuclein (Phospho-Ser129) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P37840
Expiration Date	12 months from date of receipt.