

Product Datasheet

Smad3 (Phospho-Ser423/Ser425) Rabbit mAb Antibody (orb1727815)

Description	Smad3 (Phospho-Ser423/Ser425) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P84022
Expiration Date	12 months from date of receipt.