

## Product Datasheet

# Smad5 (Phospho-Ser463/Ser465) Rabbit mAb Antibody (orb1727816)

<b>Description</b>	Smad5 (Phospho-Ser463/Ser465) Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Storage</b>	Store at -20°C. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Recombinant Rabbit
<b>Uniprot ID</b>	<a href="#">Q99717</a>
<b>Expiration Date</b>	12 months from date of receipt.