

Product Datasheet

c-Myc (Phospho-Thr58/Ser62) Rabbit mAb Antibody (orb1727827)

Description	c-Myc (Phospho-Thr58/Ser62) Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Storage	Store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.