

## Product Datasheet

### CD71 Antibody (orb1727903)

<b>Description</b>	CD71 Monoclonal Antibody [DF1513]- Biotin
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	KGI cell line
<b>Target</b>	CD71
<b>Preservatives</b>	Purified antibody (from supernatant) containing PBS + 0.09% sodium azide, conjugated to Biotin
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P02786</b>
<b>Entrez</b>	<b>7037</b>
<b>Expiration Date</b>	12 months from date of receipt.