

## Product Datasheet

### Mouse GLP-1 Ready-To-Use IHC Kit (orb1727914)

|                        |   |
|------------------------|---|
| <b>Description</b>     | orb1727914 is a ready-to-use IHC kit for staining of GLP-1. The kit provides all reagents, from anti...                                       |
| <b>Reactivity</b>      | Mouse   |
| <b>Storage</b>         | Please store components at the temperatures indicated on the individual tube labels. The kit is stable for 6 months from the date of receipt. |
| <b>Note</b>            | For research use only   |
| <b>Sample Types</b>    | PFPE tissue   |
| <b>Expiration Date</b> | Please enquire.   |