

Product Datasheet

Recombinant mouse P-Selectin, C-His (HEK293) (orb1727918)

Description	Recombinant mouse P-Selectin, C-His (HEK293)
Endotoxins	Not analyzed
Target	SELP
Preservatives	PBS (pH=7.4)
Form/Appearance	Lyophilized or Liquid
Concentration	□0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	42-709/768
Purity	> 90% as determined by SDS-PAGE
MW	110 kDa
Uniprot ID	Q01102
Expiration Date	6 months from date of receipt.