

## Product Datasheet

### Recombinant mouse TLR4 protein (orb1727920)

<b>Description</b>	Recombinant mouse TLR4 protein
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	TLR4
<b>Preservatives</b>	8M Ures(pH=8.0) with 20 mM Tris
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	> 1mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Protein Sequence</b>	26-629/835
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	70 kDa
<b>Uniprot ID</b>	<a href="#">Q9QUK6</a>
<b>Expiration Date</b>	6 months from date of receipt.