

## Product Datasheet

# Recombinant human CD34 protein-C-His (HEK293) (orb1727924)

|                         |  |
|-------------------------|--|
| <b>Description</b>      | Recombinant human CD34 protein-C-His (HEK293)                |
| <b>Endotoxins</b>       | Not analyzed   |
| <b>Target</b>           | CD34   |
| <b>Preservatives</b>    | PBS (pH=7.4)   |
| <b>Form/Appearance</b>  | Lyophilized or Liquid  |
| <b>Concentration</b>    | > 0.5mg/ml   |
| <b>Storage</b>          | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles. |
| <b>Tag</b>              | C-His  |
| <b>Note</b>             | For research use only  |
| <b>Protein Sequence</b> | 277-400/400  |
| <b>Purity</b>           | > 90% as determined by SDS-PAGE                              |
| <b>MW</b>               | 14.7 kDa   |
| <b>Uniprot ID</b>       | <b>P28906</b>  |
| <b>Expiration Date</b>  | 6 months from date of receipt.                               |