

## Product Datasheet

# Recombinant mouse CD62p protein,C-His (HEK293) (orb1727925)

<b>Description</b>	Recombinant mouse CD62p protein,C-His (HEK293)
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	SELP
<b>Preservatives</b>	PBS (pH=7.4)
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	> 0.5mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Protein Sequence</b>	455-561/768
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	13.6 kDa
<b>Uniprot ID</b>	<a href="#">Q01102</a>
<b>Expiration Date</b>	6 months from date of receipt.