

Product Datasheet

PRRSV Antibody (orb1727960)

Description	PRRSV Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant PRRSV N protein (1-128/128 aa)
Target	PRRSV
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13.8 kDa
Uniprot ID	Q6YP02
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.