


Product Datasheet

Mouse Layilin ELISA Kit (orb1728075)

Description	Mouse Layilin ELISA Kit (96 Tests). Quantitate Mouse Layilin in cell culture supernatants, serum and...
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K18-E229
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q8C351
Sensitivity	75 pg/ml
Expiration Date	Please enquire.



Mouse
Layilin
PicoKine
ELISA Kit
standar...