

Product Datasheet

Anti Mouse IL 1 β Antibody (orb1729603)

Description	Rabbit polyclonal antibody raised against recombinant truncated mouse IL1 β protein
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant truncated mouse IL1 β protein
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	NA
Hazard Information	For Research use only
Dilution Range	WB- 1:1000, ELISA 1:5000
Expiration Date	12 months from date of receipt.