

Product Datasheet

Anti Streptavidin Ab (orb1729607)

Description	Rabbit polyclonal antibody raised against recombinant full length Streptavidin protein
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length Streptavidin protein from Streptomyces avidinii
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	NA
Hazard Information	For Research use only
Dilution Range	WB-1:2000 ELISA- 1:5000
Expiration Date	12 months from date of receipt.