

Product Datasheet

Strep-tagII-Tag Antibodies (orb1729635)

Description	Rabbit polyclonal antibody to Strep-tag®II-Ta
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide Strep-tagII epitope tag, sequence: ASWSHPQFEKGA
Form/Appearance	Lyophilized
Storage	For reconstitution add 50 µl, of sterile water. Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.
Note	For research use only
Clonality	Polyclonal
Purity	Immunogen affinity purified serum, in PBS pH 7.4.
Dilution Range	1 : 1000 (WB)
Expiration Date	12 months from date of receipt.