

## Product Datasheet

### DYKDDDDK Antibody (orb1729636)

<b>Description</b>	Rabbit polyclonal antibody to DYKDDDDK (binds to Sigma FLAG®)
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide: DYKDDDDK (Sigma FLAG)
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	For reconstitution add 50 µl, of sterile water. Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity purified serum, in PBS pH 7.4.
<b>MW</b>	Depends upon fusion protein partner
<b>Uniprot ID</b>	UniProt
<b>Dilution Range</b>	1 : 1000 - 1: 5000 (WB)
<b>Expiration Date</b>	12 months from date of receipt.