

Product Datasheet

RhoA/B/C Antibody (orb1729805)

Description	RhoA,RhoB,RhoC is a small G protein of the Rho family. Regulates a signal transduction pathway linki...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Immunogen	Recombinant protein of human Rho C
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Uniprot ID	P62745 , P61586 , P08134
Entrez	387 , 389 , 388
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.