

## Product Datasheet

### PSMD14 Antibody (orb1729807)

<b>Description</b>	Metalloprotease component of the 26S proteasome that specifically cleaves 'Lys-63'-linked polyubiqui...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Immunogen</b>	A synthetic peptide of human PSMD14
<b>Preservatives</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 35 kDa; Observed MW: 35 kDa
<b>Uniprot ID</b>	<a href="#">O00487</a>
<b>Entrez</b>	<a href="#">10213</a>
<b>Dilution Range</b>	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200
<b>Expiration Date</b>	12 months from date of receipt.