

Product Datasheet

p21 Antibody (orb1729815)

Description	The tumor suppressor protein p21 Waf1/Cip1 acts as an inhibitor of cell cycle progression. It functi...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	A synthetic peptide of human p21
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 18 kDa; Observed MW: 21 kDa
Uniprot ID	P38936
Entrez	1026
Dilution Range	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.