

## Product Datasheet

### RhoA/B/C Antibody (orb1729857)

<b>Description</b>	RhoA,RhoB,RhoC is a small G protein of the Rho family. Regulates a signal transduction pathway linki...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IHC, WB
<b>Immunogen</b>	A synthesized peptide derived from human Rho A + B + C
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 22 kDa; Observed MW: 22 kDa
<b>Uniprot ID</b>	<a href="#">P62745</a> , <a href="#">P61586</a> , <a href="#">P08134</a>
<b>Entrez</b>	<a href="#">387</a> , <a href="#">389</a> , <a href="#">388</a>
<b>Dilution Range</b>	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
<b>Expiration Date</b>	12 months from date of receipt.