

Product Datasheet

BTN3A Antibody (orb1729863)

Description	The butyrophilin (BTN) genes are a group of major histocompatibility complex (MHC)-associated genes ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide from human protein at AA range: 41-90
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P78410 , O00481 , O00478
Entrez	10384 , 11118 , 11119
Dilution Range	IHC: 1/50-1/100; ELISA: 1/10000
Expiration Date	12 months from date of receipt.