

## Product Datasheet

### RPA3 Antibody (orb1729892)

<b>Description</b>	As part of the heterotrimeric replication protein A complex (RPA/RP-A), binds and stabilizes single-...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide from human protein at AA range: 30-80
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 14 kDa; Observed MW: 14 kDa
<b>Uniprot ID</b>	<a href="#">P35244</a>
<b>Entrez</b>	<a href="#">6119</a>
<b>Dilution Range</b>	WB: 1/500-1/1000; ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.