

## Product Datasheet

### FGF18 Antibody (orb1730025)

<b>Description</b>	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF famil...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FGF18. AA range:158-207
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 24 kDa; Observed MW: 28 kDa
<b>Uniprot ID</b>	<a href="#">076093</a>
<b>Entrez</b>	<a href="#">8817</a>
<b>Dilution Range</b>	WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.