

---

## Product Datasheet

### Ubiquitin Conjugating Enzyme E2 A/B Antibody (orb1730059)

**Description**

Accepts ubiquitin from the E1

attachment to other proteins. In as...

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IHC, WB

**Immunogen**

Synthesized peptide derived from the C-terminal region of human UBE2A/B.

**Preservatives**

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

Rabbit IgG

**Clonality**

Polyclonal

**MW**

Calculated MW: 17 kDa; Observed MW: 17 kDa

**Uniprot ID**[P49459](#), [P63146](#)**Entrez**[7319](#), [7320](#)**Dilution Range**

WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000

**Expiration Date**

12 months from date of receipt.