

Product Datasheet

p21 Antibody (orb1730076)

Description	The tumor suppressor protein p21 Waf1/Cip1 acts as an inhibitor of cell cycle progression. It functi...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human p21 Cip1. AA range:111-160
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 18 kDa; Observed MW: 26 kDa
Uniprot ID	P38936
Entrez	1026
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000
Expiration Date	12 months from date of receipt.