

Product Datasheet

DLC1 Antibody (orb1730087)

Description	Functions as a GTPase-activating protein for the small GTPases RHOA, RHOB, RHOC and CDC42, terminati...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide of human DLC1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 171 kDa; Observed MW: 171,127,123,114 kDa
Uniprot ID	Q96QB1
Entrez	10395
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.