

Product Datasheet

NMDAR2B Antibody (orb1730088)

Description	Component of NMDA receptor complexes that function as heterotetrameric, ligand-gated ion channels wi...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	Synthetic peptide of human GRIN2B
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q13224
Entrez	2904
Dilution Range	ICC: 1/100-1/200
Expiration Date	12 months from date of receipt.