

Product Datasheet

HDAC5 Antibody (orb1730142)

Description	Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human HDAC5. AA range:464-513
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Uniprot ID	Q9UQL6
Entrez	10014
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; ELISA: 1/10000
Expiration Date	12 months from date of receipt.