

## Product Datasheet

### Fragilis Antibody (orb1730303)

<b>Description</b>	IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human p...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IHC, WB
<b>Immunogen</b>	A synthesized peptide derived from human Fragilis
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 15 kDa; Observed MW: 15 kDa
<b>Uniprot ID</b>	<a href="#">Q01628</a>
<b>Entrez</b>	<a href="#">10410</a>
<b>Dilution Range</b>	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
<b>Expiration Date</b>	12 months from date of receipt.