

Product Datasheet

HDAC8 Antibody (orb1730430)

Description	Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	A synthesized peptide derived from human HDAC8
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	Q9BY41
Entrez	55869
Dilution Range	IP: 1/20; WB: 1/500-1/1000; RCM: 1/50-1/100
Expiration Date	12 months from date of receipt.